

Data sheet

Ultra-Low Freezer ZLN-UT 200



The photo above is for reference only, may show additional options not included in standard equipment.
The real appearance, particularly color and structure of the material may differ from the ones presented in the photo.


TECHNICAL DATA	
air convection	natural
chamber capacity [l]	237
working capacity [l]	237
controller	microprocessor with external LCD graphic display
TEMPERATURE	
temperature range [°C]	-86...-40
temperature resolution every ... [°C]	0,1
over temperature protection	-
CHAMBER	
door type	solid
interior	
COMF	stainless steel to DIN 1.4016
PREM	stainless steel to DIN 1.4301
housing	
COMF	powder coated sheet
PREM	powder coated sheet

overall dims [mm] /1/

width A	850
height B	1620
depth C	950

internal dims [mm]

width D	520
height E	830
depth F	550
depth G	240



number of internal chambers	3
max shelf workload [kg]	10
max unit workload [kg]	65
weight [kg]	180

ELECTRICAL PARAMETERS

voltage*	230V 50Hz
nominal power [W]	670
warranty	24 months
manufacturer	POL-EKO-APARATURA

all the above technical data refer to standard units (without optional accessories)

* - 230V 60Hz, 115V 60Hz also available

1 - depth doesn't include 50 mm of power cable

OPTIONS AND ACCESSORIES

