

## Data sheet

### Ultra-Low Freezer ZLN-UT 300



The photo above is for reference only, may show additional options not included in standard equipment.  
The real appearance, particularly color and structure of the material may differ from the ones presented in the photo.

#### TECHNICAL DATA

air convection	natural
chamber capacity [l]	326
working capacity [l]	326
controller	microprocessor with external LCD graphic display

#### TEMPERATURE

temperature range [°C]	-86...-40
temperature resolution every ... [°C]	0,1
over temperature protection	-

#### CHAMBER

door type	solid
-----------	-------

#### interior

COMF	stainless steel to DIN 1.4016
PREM	stainless steel to DIN 1.4301

#### housing


COMF	powder coated sheet
PREM	powder coated sheet

**overall dims [mm] /1/**

width A	850
height B	1910
depth C	950

**internal dims [mm]**

width D	520
height E	1140
depth F	550
depth G	240



number of internal chambers	4
max shelf workload [kg]	10
max unit workload [kg]	80
weight [kg]	200

**ELECTRICAL PARAMETERS**

voltage*	230V 50Hz
nominal power [W]	670
warranty	24 months
manufacturer	POL-EKO-APARATURA

all the above technical data refer to standard units (without optional accessories)

\* - 230V 60Hz, 115V 60Hz also available

1 - depth doesn't include 50 mm of power cable

**OPTIONS AND ACCESSORIES**

