

TROLLEYS™

- Entirely made of Stainless Steel (AISI 304)
- Ø 125 mm casters,
two of which can be equipped with brakes
(on demand)
- Available in various models



LOGISTICS



TECNIPLAST
innovation through passion™

1/Laboratory Table

Suitable for various needs: from cage changing to animal weighing. Supplied with two shelves with rounded edges for complete safety.

2/Cage transport trolley

Suitable for the transport of stacked plastic cages, bedding and feed bags, etc.

Ergonomic shape ensures easy, safe movement. Supplied with a Stainless Steel (AISI 304) shelf and front fenders.

3/Feed container

The trolley has two separate 80 lt. bins to take two different types of feed. Separate lids seal off the two bins avoiding contamination from the outside. The inner floor of the two bins is made of perforated Stainless Steel sheet to minimise build up of dust during normal use. The dust residue falls into removable collection trays positioned below.

[The following models are also available: FEED CONTAINERS with only one compartment - GLASSWARE TRANSPORT UNITS - OTHER TROLLEYS on demand].



LOGISTICS

Trolleys

Laboratory table **1**



Cage transport trolley **2**



Feed container **3**



TECHNICAL SHEET

Laboratory Table

Cod. ACCACG01 Working space - mm: 950 x 522
Overall dimensions (L x W x H) - mm: 1100 x 565 x 870

Cod. ACCACG02 Working space - mm: 1200 x 522
Overall dimensions (L x W x H) - mm: 1350 x 565 x 870

Cage Transport Trolley

Cod. ACCAPA02 Overall dimensions (L x W x H) - mm: 970 x 650 x 182

Feed Container

Cod. ACCAPM00 Overall dimensions (L x W x H) - mm: 1100 x 565 x 870
Cod. ACCAPM02 Overall dimensions (L x W x H) - mm: 775 x 565 x 870 (1 bin)

Bedding Container

Similar to the feed container without perforated inner floor (160 lt. capacity)

Cod. ACCAPS00 Overall dimensions (L x W x H) - mm: 1100 x 565 x 870 (2 bins)
Cod. ACCAPS01 Overall dimensions (L x W x H) - mm: 1100 x 565 x 870 (1 bin)

TECNIPLAST

HEADQUARTERS

ITALY | Tel. +39 0332 809711 • www.tecniplast.it • E-mail: tecnicom@tecniplast.it

SISTER COMPANIES

AUSTRALIA/NEW ZEALAND | Tel. + 61 2 8845 6500 • www.tecniplast.it • E-mail: info@tecniplast.com.au

FRANCE | Tel. 04 72 52 94 41 • www.tecniplast.fr • E-mail: info@tecniplast.fr

GERMANY | Tel. 08805 921320 • www.tecniplast.de • E-mail: info@tecniplast.de

UNITED KINGDOM | Tel. 0845 0504556 • www.tecniplast.it • E-mail: info@tecniplastuk.com

JAPAN | Tel. + 03 5770 5375 • www.tecniplast.it • E-mail: info@tecniplastjapan.co.jp

USA/CANADA | Toll Free: 877.669.2243 • www.tecniplastusa.com • E-mail: info@tecniplastusa.com

To find your local distributor, please visit www.tecniplast.it



ISO-9001 - Cert. n° 0875



ISO-14001 - Cert. n° 0256A/0



Tecniplast reserves the right to change or modify product and or specifications without notice or obligation. Trolleys is a Tecniplast trademark.

CAT015/ING - Rev. March 2011