

- The size you have asked for is here!
- A 13-liter water capacity spawning system that allows you to collect thousands of embryos based on your research needs
- Fits on a benchtop. Easy to handle



## iSPAWN - S

Because you've asked for it!

Designed to efficiently collect a large number of synchronized Zebrafish embryos in brief and defined windows of time, with minimal labor and space



**AQUATIC SOLUTIONS**

• **Synchronized embryo collection that suits your needs**

iSPAWN-S is the new 13-liter capacity spawning system that allows researchers to breed fish and collect embryos in batches of thousands at time.

• **Natural spawning behaviour**

iSPAWN-S capitalizes on the fish natural tendency to spawn in shallow water. The innovative design promotes mating behaviour that results in intense bursts of embryo production.

• **Easy to use**

iSPAWN-S was designed with convenience and flexibility in mind. It fits on a benchtop and is easy to handle. During the breeding protocol, males and females are sequentially added to the system in the evening and kept separated by a divider. The next day, when the embryos are required, the divider is removed and the spawning platform is lifted to initiate the breeding event.

• **Scientifically proven**

iSPAWN-S has been successfully tested by experts in the most active zebrafish labs in the world.

• **Fish welfare**

iSPAWN-S is completely transparent for maximum visibility to facilitate inspection of the fish.

• **Flexible solution**

iSPAWN-S offers the flexibility to produce embryos based on your lab needs by allowing single or multiple configurations, independently from your research focus.

**By introducing iSPAWN-S, we offer the market widest selection of breeding solutions, to fit all your research needs**

For information about iSPAWN-S and other breeding systems contact your local representative or visit [www.aquaticsolutions.it](http://www.aquaticsolutions.it)



**Breeding tank**  
0.7-liter



**Breeding tank**  
1-liter



**Breeding tank**  
2-liter



**iSPAWN-S**  
13-liter



**iSPAWN**  
60-liter



**TECNIPLAST**

HEADQUARTERS  
ITALY | Tel. +39 0332 809711 • [www.tecniplast.it](http://www.tecniplast.it) • E-mail: [tecnicom@tecniplast.it](mailto:tecnicom@tecniplast.it)

SISTER COMPANIES

AUSTRALIA/NEW ZEALAND | Tel. + 61 2 8845 6500 • [www.tecniplast.it](http://www.tecniplast.it) • E-mail: [info@tecniplast.com.au](mailto:info@tecniplast.com.au)

CHINA | Tel. + 86 (0)21 50810920 • [www.tecniplast.cn](http://www.tecniplast.cn) • E-mail: [tecniplastchina@tecniplast.it](mailto:tecniplastchina@tecniplast.it)

FRANCE | Tel. +33 (0)4 72 52 94 41 • [www.tecniplast.fr](http://www.tecniplast.fr) • E-mail: [info@tecniplast.fr](mailto:info@tecniplast.fr)

GERMANY | Tel. +49 (0)8805 921320 • [www.tecniplast.de](http://www.tecniplast.de) • E-mail: [info@tecniplast.de](mailto:info@tecniplast.de)

UNITED KINGDOM | Tel. + 44 (0)845 0504556 • [www.tecniplastuk.com](http://www.tecniplastuk.com) • E-mail: [info@tecniplastuk.com](mailto:info@tecniplastuk.com)

JAPAN | Tel. +81 (0)3 5770 5375 • [www.tecniplast.jp](http://www.tecniplast.jp) • E-mail: [info@tecniplastjapan.co.jp](mailto:info@tecniplastjapan.co.jp)

USA | Toll Free: +1 877.669.2243 • [www.tecniplastusa.com](http://www.tecniplastusa.com) • E-mail: [info@tecniplastusa.com](mailto:info@tecniplastusa.com)

CANADA | Toll Free: +1 855.347.8718 • [www.tecniplastcanada.com](http://www.tecniplastcanada.com) • E-mail: [info@tecniplastcanada.com](mailto:info@tecniplastcanada.com)

To find your local distributor, please visit [www.tecniplast.it](http://www.tecniplast.it)

Tecniplast reserves the right to change or modify product and/or specifications without notice or obligation.

CAT 154/ING - Rev. November 2015